



THE UNIVERSITY OF
CHICAGO
BIOLOGICAL
SCIENCES

Flow Basics 2.1: The Basic Staining Protocol and How to Improve Reproducibility

May 2020


Laura K Johnston, PhD

Scientific Associate Director

Cytometry and Antibody Technology Facility

University of Chicago

The Flow Basics 2.0 Series



Understanding Flow Cytometry Experiments to Get Better Results

- For all scientific experiments the best data is achieved by optimization and consistency!
- This course will go over:
 - A basic staining protocol for surface markers
 - How to improve consistency in sample staining

Staining Protocol

Staining Tubes and Buffers

Staining tube:


- Eppendorf, 5 mL FACS tube, or 96-well round bottom plate
- Choice depends on preference. Main difference is wash volume:



Buffer


- Choice depends on your cell type
- FACS buffer options:
 - 1X PBS (Ca/Mg²⁺ free, no phenol red!)
 - Either FBS (1-10%) or BSA (0.1-1%)
 - Commonly used: 1% FBS in 1X PBS (consider filtering FBS to remove debris, especially for cell sorting)
 - Other additives: EDTA (0.5-5mM) if you have very sticky cells, sodium azide

Basic staining protocol




Blocking to reduce nonspecific and Fc-mediated staining

- Immune cells have Fc receptors that bind IgG antibodies. Most antibodies used for flow cytometry are IgG
 - Fc Block inhibits IgG antibodies from binding to Fc receptors
 - In mice, Fc block is anti-CD16/CD32
 - Serum could also be used
 - E.g. Rabbit serum for rabbit antibodies, goat serum for goat antibodies
 - Also BSA/FBS in staining buffer helps



Basic staining protocol

Staining Surface Markers with Fluorescent Antibodies



Antibody Staining is Affected by Five Factors

1. Time
2. Temperature
3. Volume
4. Cell number
5. Antibody amount

Antibody Staining is Affected by Five Factors

1. Time

2. Temperature

3. Volume

4. Cell number

5. Antibody amount

Time and temperature

- Typically 30 minutes at 4°C
- Time is often between 15-60 minutes, sometimes overnight
- Temperature is often 4°C or room temp
- Some antibodies/dyes require something else – read manufacturer protocols!

Pick a time and temperature and don't vary the protocol between experiments!

Antibody Staining is Affected by Five Factors

1. Time
2. Temperature
- 3. Volume**
4. Cell number
5. Antibody amount

Staining Volume

- Typically 100 μL
- Some antibodies/dyes require a different volume
- Flow Basics 2.2 will discuss when changing the staining volume should be considered

Antibody Staining is Affected by Five Factors

1. Time
2. Temperature
3. Volume
4. Cell number
5. Antibody amount

Cell Number and Antibody Amount

- Typical starting point is 1×10^6 cells and 0.1 μg of antibody (usually 0.5-1 μL)
- Will discuss further in Flow Basics 2.2

Pipetting – another source of variability


- Avoid pipetting small volumes of antibody
 - If you must pipette less than 1-2 μL , make a 1:10 dilution (pipette 5 μL instead of 0.5 μL)
 - Make a master mix of your antibody cocktail to ensure all of your samples are stained with the same amount of antibody
- Simplify your pipetting
 - You can add Fc block to your cells in 50 μL
 - Add FACS buffer to master mix to bring up to 50 μL
 - To stain, add 50 μL cells and 50 μL master mix for a total staining volume of 100 μL

	1X	10X
FITC	1	10
APC	5	50
PE	0.5	5
FACS Buffer	43.5	435
TOTAL	50	500

Always protect fluorophores from light!

- Once fluorescent antibodies have been added to the cells, be sure to keep the cells in the dark as much as possible
- Most common: Cover tubes with foil and incubate in a dark place (fridge for 4°C or a drawer or cabinet for room temperature)
- Note: cells with fluorescent protein(s) should always be protected from light

Basic staining protocol



Wash and centrifuge

- Volume depends on tube – 1 mL, 1-3 mL, or 0.2 mL
- Wash buffer with protein (BSA/FBS) helps prevent cell loss compared to PBS without protein
- Aspirate with dumping or pipetting, careful with vacuum

Resuspending Cells to Run on a Cytometer

- If you are new to flow cytometry, it is safer to resuspend your cells (especially controls) in a larger volume. This will allow you to properly set the instrument settings without running out of sample
 - As you become more comfortable with the cytometer, you can use more concentrated controls so that they run faster on the instrument
- The volume to resuspend cells in depends on the instrument
 - Fortessas
 - 150 μL **absolute minimum**, 300-500 μL minimum for new users, 150-300 μL for advanced users
 - In general, run samples at less than 35,000 events/sec
 - Attune: see chart [here](#)
 - Sorters
 - 500 μL **absolute minimum**
 - 70um tip (lymphocytes, splenocytes, small cells): 20-25 million cells per mL
 - 100um tip (most cell culture, lung and liver cells, etc): 10-15million cells per mL


Basic Staining Protocol: Additional Notes

- Use a viability stain to remove dead cells from analysis
- Using primary antibodies vs primary and secondary antibodies?
- Should a fixative be used?

Viability Dyes


- Dead cells will:
 - Bind some antibodies non-specifically
 - Skew statistics
- DNA-binding (DAPI, PI, etc.)
 - Non-fixable
 - Stain **after antibodies**, right before acquisition
- Amine-binding
 - Fixable (Biolegend Zombie dyes or ThermoFisher Fixable Live/Dead dyes)
 - Stain **before antibodies** and **before** fixation (fixation not required)
 - Staining buffer must be free of protein – read the manufacturers protocol!!

Zombie Dyes



*DRAQ7™ is a trademark of Biostatus Limited.

Helix NP™ NIR, DRAQ7™, PI and 7-AAD



Product Links: [Helix NP™ NIR](#), [DRAQ7™](#), [PI](#), [7-AAD](#)

More info: <https://www.biolegend.com/en-us/live-dead>

Using Primary and Secondary Antibodies

- Flow cytometrists prefer to use fluorophore-conjugated primary antibodies as opposed to unconjugated primary antibodies and fluorophore-conjugated secondary antibodies
 - Using a primary antibody instead of a primary and secondary:
 - Only requires one incubation time and wash steps instead of two
 - Requires less optimization (one antibody instead of two)
 - Makes the staining protocol easier
 - Avoids chances of nonspecific/background staining from the secondary antibody
- You are welcome to use primary and secondary antibodies if you want – just titrate both antibodies!

Fixatives

- Fresh cells (BSL1) can be run on all of our instruments
- Cells with infectious particles (BSL2 and above) are required to be fixed prior to running on the benchtop analyzers
- There are many options for fixatives: 1-4% paraformaldehyde, methanol, etc.
- Concerns with fixation:
 - Fixing cells before antibody staining may alter the cell's epitope, leading to altered or no staining
 - Fixing cells after antibody staining may decrease the signal intensity of some fluorophores
 - Fixing cells could increase cell autofluorescence, which may decrease resolution of markers

Resources for Fixation

- <https://www.biolegend.com/en-us/blog/fix-now-fix-later-considerations-for-the-use-of-paraformaldehyde-fixation-in-flow-cytometry>
- <https://bitesizebio.com/22141/fixation-and-flow-cytometry/>

Stay Tuned for the Rest of the Flow Basics 2.0 Series

