


Instrument Configuration



Instrument: NovoCyte Penteon

Optical Configuration: U7V7B6Y6R4



Par...	Excit...	Detectio...	Default Alias
B525	488	525/45	FITC
B586	488	586/20	EYFP
B615	488	615/20	PI
B667	488	667/30	PerCP
B695	488	695/40	PerCP-Cy5.5
B725	488	725/40	PerCP-eFluor 7...
R667	637	667/30	APC
R695	637	695/40	Alexa Fluor 680
R725	637	725/40	Alexa Fluor 700
R780	637	780/60	APC-Cy7
UV445	349	445/45	BUV395
UV525	349	525/45	BUV496
UV586	349	586/20	BUV563
UV615	349	615/20	BUV615
UV667	349	667/30	BUV661
UV725	349	725/40	BUV737
UV780	349	780/60	BUV805
V445	405	445/45	Pacific Blue
V525	405	525/45	AmCyan
V586	405	586/20	Pacific Orange
V615	405	615/20	Qdot 605
V667	405	667/30	Qdot 655
V725	405	725/40	BV711
V780	405	780/60	Qdot 800

NovoSampler Q installed

OK

Cancel